

Debulking Buffer

Model # R2137

WAVESENSE

Description:

This product is a pH 7.5 buffered non-ionic detergent solution, intended for In Vitro Diagnostic Use.

Supplied As:

Catalog #	Contains
R2137-1	100 mL

A proprietary buffered solution containing surfactant and preservatives.

Storage:

This product is stable when stored at 18 - 25°C. DO NOT FREEZE. Refer to product label for expiration date.

Material Safety Data:

When handling this material, Standard Laboratory Practices should be followed. This material's chemical, physical and toxicological properties have not been thoroughly investigated. Contains Sodium Azide as a preservative. Although, the quantity of Sodium Azide (<0.9%) is very small, measures should be taken to avoid skin and eye contact, inhalation and ingestion. Sodium Azide (NaN₃) may react with lead copper plumbing to form potential explosive metal oxides. Upon disposal, flush with a large volume of water to prevent Azide build-up.

Product Specification Sheet



15339 Barranca Pkwy
Irvine, CA 92618 USA
www.WaveSense.net

Toll Free: 800.807.7760
Phone: 949.341.1980
Fax: 949.341.1982
Contact@WaveSense.net

ML100303-22
© 02/08/13

Page 1 of 1